


Secura-sSAY



on a scale of 5:1

Technical Information

Type	Posterior Chamber IOL with Blue Light Protection (Three-Piece, Foldable) for Sulcus Fixation
Optic Material	Silicone Elastomer with UV-Absorber and Yellow Colouring Agent
Haptic Material	High Molecular PMMA
Haptic Shape	Modified C-Loop
Optic Shape	Spheric Anterior Surface, Aspheric Posterior Surface (Aberration Correcting)*
	Biconvex
Diopter Range	7.0 to 31.0 D
Diopter Steps	0.5 D
A-Constant (est.)**	118.0

HumanOptics AG

Spardorferstraße 150
91054 Erlangen

Phone: +49 (0) 9131 / 50665-0
Fax: +49 (0) 9131 / 50665-90

www.humanoptics.com
mail@humanoptics.com

03/2011

* The word „aberration“ as used in this document is interpreted as „spherical aberration“

** For current status and more details please refer to www.humanoptics.com